
SpermoLyzer (Computer Assisted Semen Analysis)

The system follows WHO (2010) strict criteria for motility patterns & morphometric assessment of human semen.

Name	محمد عبدالله محمد الطاينس	Case no.	202210245	
Age	25	Date	24/10/2022	
Lab ID	110737	Ref. by	بلقیس شریان∖.Dr	

	Physical Properties		
Abstinence	4	2 - 5 days	
Collection	Masturbation at lab		
Volume	6 >= 1.5 ml		
рН	7.8	>= 7.2	
Color	Grey Opalescent		
Odor	Normal		
Viscosity	High		
Liquefaction time	110 15 - 60		
Liquefaction State	Complete		



Test Result				
Test	Result	Status	Reference values	
Concentration	18.4	Passed	>= 15 million/ml	
Total sperm count	110.4	Passed	>= 39 million/ejaculate	
Progressive motility (PR)	47.27	Passed	>= 32 %	
Rapid PR	29.04			
Slow PR	18.23			
Total motility (PR+NP)	70.94	Passed	>=40 %	
Morphology index	9	Passed	>= 4 %	

		Agglutination grade
Agglutination	Nil	
Agglutination grade	-	
Non-specific Aggregation	Nil	
	Other Cells	
White blood cells	1 - 2	< 5 /H.P.F
Red blood cells		< 5 /H.P.F
Spermatogenic cells	0 - 1	/HPF
Epithelial cells		/HPF

Normal seminal profile

The semen passed the WHO strict criteria for count and motility.

SpermoLyzer (Computer Assisted Semen Analysis) Motility report

Name	محمد عبدالله محمد الطايفي	Lab ID	110737	

Dynamic parameters (I)					
Classification Conc. (million/ml) Total number (million) Percentage (%)					
Tested sperms	18.4	110.4	100		
Total motility (PR+NP)	13.05	78.3	70.94		
Progressive motility (PR)	8.7	52.2	47.27		
Non progressive (NP)	4.36	26.16	23.67		
Immotile (IM)	5.35	32.1	29.06		

Progressive motility (PR) Spermatozoa moving actively either linearly or in a large circle regardless of speed

Non progressive (NP) All other patterns of motility with absence of progression, i.e. swimming in small circle, the flagellar force hardly displacing the head, or when only a flagellar beat can be observed.

Immotile (IM) No movement

Dynamic parameters (II)					
VCL (μm/s) 28.64 VSL (μm/s) 13.32 VAP (μm/s) 17.63					
MAD (°)	43.36	ALH (µm)	5.02	BCF (Hz)	4.13
LIN (%)	46.18	WOB (%)	61.4	STR (%)	75.12

VCL Curvilinear velocity	VSL Straight line velocity	VAP Average path velocity
MAD Mean angular degree	ALH Amplitude of lateral head displacement	BCF Beat-cross frequency
LIN Linearity (VSL/VCL)	WOB Wobble (VAP/VCL)	STR Straightness (VSL/VAP)




Diagram of sperm tracks and velocities

Legend

SpermoLyzer (Computer Assisted Semen Analysis) Charts Report

